

Products List

Fluorinated Alcohols

2N004	2,2,2-Trifluoroethanol [75-89-8] M.W. 100.0 b.p. 73.6°C d ²⁰ 1.384	CF ₃ CH ₂ OH (TFEA)
4N032	4,4,4-Trifluoro-1-butanol [461-18-7] M.W. 128.1 b.p. 125°C d ²⁰ 1.220	CF ₃ CH ₂ CH ₂ CH ₂ OH
4N033	2-(2,2,2-Trifluoroethoxy)ethanol [2358-54-5] M.W. 144.1 b.p. 76°C/50mmHg d ²⁴ 1.29	CF ₃ CH ₂ OCH ₂ CH ₂ OH
3N026	2,2,2-Trifluoroethyl methanesulfonate [25236-64-0] M.W. 178.1 b.p. 97-99°C/35mmHg d ²⁰ 1.513	CH ₃ SO ₂ OCH ₂ CF ₃

Fluorinated Olefins

3E011	2-Bromo-3,3,3-trifluoropropene [1514-82-5] M.W. 174.9 b.p. 35°C d ¹⁵ 1.686	CF ₃ CBr=CH ₂
4E004	4-Bromo-3,3,4,4-tetrafluoro-1-butene [18599-22-9] M.W. 207.0 b.p. 63°C d ²⁰ 1.622	CF ₂ BrCF ₂ CH=CH ₂

Iodofluorocarbons (IFCs)

1N007	Iodotrifluoromethane [2314-97-8] M.W. 195.9 b.p. -22.5°C d ²⁰ 2.316	CF ₃ I
2N019	Iodopentafluoroethane [354-64-3] M.W. 245.9 b.p. 13°C d ²⁸ 2.072	CF ₃ CF ₂ I
2N036	2-Iodo-1,1,1-trifluoroethane [353-83-3] M.W. 209.9 b.p. 54.8°C d ²⁰ 2.130	CF ₃ CH ₂ I
4N031	4-Iodo-1,1,1-trifluorobutane [461-17-6] M.W. 238.0 b.p. 128°C d ²⁰ 1.856	CF ₃ CH ₂ CH ₂ CH ₂ I

Trifluoromethylating Agent

4N024	(Trifluoromethyl)trimethylsilane [81290-20-2] M.W. 142.2 b.p. 57°C d ²⁰ 0.9626	CF ₃ Si(CH ₃) ₃
-------	---	---










TOSOH F-TECH, INC.





TOSOH

Products List

Fluorinated carboxylic Acids

5E011	2,2,2-Trifluoroethyl acrylate [407-47-6] M.W. 154.1 b.p. 46°C/125mmHg d ²⁵ 1.216		
6E009	2,2,2-Trifluoroethyl methacrylate [352-87-4] M.W. 168.1 b.p. 107°C d ²⁰ 1.180		(Fluorester)
4E005	α-(Trifluoromethyl)acrylic acid [381-98-6] M.W. 140.1 b.p. 86°C/35mmHg m.p. 53°C		
8E019	tert-Butyl α-(Trifluoromethyl)acrylate [105935-24-8] M.W. 196.2 b.p. 37°C/23mmHg d ²⁰ 1.097		
6E030	Trifluoroethyl α-(trifluoromethyl)acrylate [91520-39-7] M.W. 222.1 b.p. 69°C/121mmHg d ²⁰ 1.435		
6E031	Ethyl 4,4,4-trifluorocrotonate [25597-16-4] M.W. 168.1 b.p. 114-115°C d ²⁵ 1.125		
6N035	Ethyl 2-fluoroacetoacetate [1522-41-4] M.W. 148.1 b.p. 183°C d ²⁵ 1.181		

Fluorinated Heterocycles

H001	5-(Trifluoromethyl)uracil [54-20-6] M.W. 180.1 m.p. 241~245°C (Dec.)	
H025	6-Trifluoromethyl-4-hydroxypyrimidine [1546-78-7] M.W. 164.1 m.p. 168~170°C	
H026	2,4-Dichloro-5-trifluoromethylpyrimidine [3932-97-6] M.W. 217.0 m.p. 61°C/8mmHg	
H027	N-tert-Butoxycarbonyl-(2S, 4S)-4-fluoroproline [203866-13-1] M.W. 253.8 m.p. 155~159°C	



TOSOH F-TECH, INC.